

### Applications of Ion-Selective Electrodes

	Ammonium	Barium	Bromide	Calcium	Cadmium	Carbonate	Chloride	Cupric	Cyanide	Fluoride	Iodide	Lead	Mercury	Nitrate	Nitrite	Perchlorate	Potassium	Silver	Sulphide	Thiocyanate
Baby Food	X			X			X			X	X			X	X					
Beer	X			X		X														
Biological Samples	X		X							X	X									
Blood										X							X			
Boiler Feedwater	X			X			X			X										
Cheese				X			X													
Drinking Water	X		X	X		X	X			X	X			X	X		X			
Environmental	X	X			X		X		X			X	X	X	X				X	
Explosives																X				
Faeces	X																X			
Feedstuffs				X							X									
Fertilizers	X			X										X			X			
Fish Bonds/Tanks	X																			
Fruits				X										X	X		X			
Gastric Juices				X			X													
Ground Water	X			X			X							X	X		X			
Liquors				X			X		X										X	
Milk				X			X			X	X						X			
Mineral Water				X		X					X									
Paper Pulp				X			X												X	
Photo Fixing Solution																		X		
Pickling Baths							X							X	X					
Plant Tissue			X				X							X			X			
Plating Baths	X	X			X		X	X	X					X				X		X
Rain Water										X										
Saliva										X										
Sea Water	X						X							X						
Serum/Plasma				X			X										X			
Sewerage Sludge	X	X						X												
Soft Drinks				X		X				X										
Soils	X		X				X			X										
Solid Propellants																				
Stack Gases	X								X											
Surface Water	X						X										X			
Sweat							X													
Teeth & Bone										X										
Toothpaste										X										
Urine	X									X							X			
Vegetables							X										X			
Waste Water	X	X	X	X			X		X								X		X	
Wine	X		X			X				X							X			